



New Psychoactive Substances - Are We Doing Enough?

July 2014

Science
for a safer world



LGC Classification: NOT PROTECTIVELY MARKED



Background

New Psychoactive Substances (NPS)



NPS (often called legal highs) cover a wide range of substances:


1. Products with names which give no indication of what they contain, Vertex Pirate Edition
2. Named and specific substances which are designed to be similar chemically and/or pharmacologically to known specific controlled drugs, for example many synthetic cannabinoids
3. Substances related to medicines, such as 2-MeO-2-deschloroketamine
4. Herbal and fungal materials or their extracts, salvia



Increasing Prevalence



- New psychoactive substances reported are increasing (EMCDDA)




■ New psychoactive substances reported



Legal Status of Legal Highs

- Legislation is slow to react to new drugs
- Many NPS are not controlled
- Fast track system – ‘Temporary Class Drugs’ in place to cover supply of most serious compounds
- Issued for benzofury (5/6 APB), NBOMe compounds and methoxetamine



Hot of the press:- Queen's Speech in UK at opening of new parliament



- New UK government approach to Legal Highs

The main elements of the Bill are:

- The Bill would make it an offence to produce, supply, offer to supply, possess with intent to supply, import or export psychoactive substances; that is, any substance intended for human consumption that is capable of producing a psychoactive effect. The maximum sentence would be seven years' imprisonment.

LGC Classification: NOT PROTECTIVELY MARKED

Risks



- Legal does **NOT** mean safe
- Most (excl prescription) not tested in humans
- Health risks (short and long term) are often not understood
- Impairment outcomes are often not understood
- Not guaranteed that you are taking what you think you are - hence possible unexpected consequences



Workplace Analytical

Traditional Techniques



- Workplace drug testing
- Screened on clinical analysers or using point of care tests (POCT) using immunoassay
- Limited panel of substances
- High throughput
- Cheap
- New reagents can take a while to develop



Alternative Approach



- High resolution accurate mass (HRAM) LCMS to cover much broader range of analytes.
- Quick and easy to update analyte coverage
- Expensive capital purchase
- Requires high level of technical expertise





Findings

City Centre Study

'Night time' economy study:

1. Pooled urine samples from urinals
2. Taken from night clubs and city centre
3. Alternative approach as opposed to self reporting
4. HRAM screening
5. 164 different substances found, including synthetic steroids, drugs of abuse, prescription and over the counter medication, and NPS



Compounds detected




1,4 methoxyphenylpiperazine	chlordiazepoxide	gabapentin	morphine	propranolol
4 methylethcathinone	chloroquine	guaifenesin	Nandrolone	pyrimethamine
4-ethylmethcathinone	chlorpheniramine	HMMA	naproxen	quetiapine
4-fluoroephedrine	Chlorthalidone	hordenine	nefopam	quinine
4-methylmethcathinone	cimetidine	hydrochlorthiazide	nevirapine	ranitidine
5/6-APB	citalopram	Ibuprofen	nicotine	risperidone
Ambroxol	clenbuterol	isometheptene	niflumic acid	salbutamol
Amisulpiride	clobazam	ketamine	nordazepam	sertaline
Amitriptyline	clomipramine	Ketoprofen	nortriptyline	sildenafil
Amphetamine	clozapine	lamotrigine	noscipine	sotalol
Ampyrone	codeine	lansoprazole	Olanzapine	Stanozolol
Arecoline	cotinine	Lidocaine	omeprazole	sulpiride
atenolol	cyclobenzaprine	MDA	opipramol	tamoxifen
Benzocaine	DEET	MDMA	oripavine	temazepam
Benzoyllecgonine	desloratidine	Meclofenamic acid	orphenadrine	tetracycline
Benzydamine	dextromethorphan	mefenamic acid	oxazepam	tetramisole
bisoprolol	dextrorphan	methadone	oxprenolol	theobromine
bromhexine	Diazepam	methamphetamine	oxycodone	theophylline
Buprenorphine	diclofenac	methcathinone	oxytetracycline	tramadol
buprenorphine-nor	dihydrocodeine	methiopropamine	papaverine	Trazodone
Caffeine	dihydromorphine	methoxetamine	paracetamol	trenbolone
Camphor	diphenhydramine	methylhexaneamine	paroxetine	trifluoromethylphenylpiperazine
Capsaicin	dipyridamole	methylphenidate	pentedrone	trimethoprim
Carbamazepine	doxycycline	metoclopramid	pheniramine	vardenafil
Carboxy THC	enalapril	metoprolol	Phenylpropanolamine	venlafaxine
Cathine	ephedrine/pseudoephedrine	metronidazole	Pholcodine	xylometazoline
Cathinone	ethylphenidate	midazolam	pipradrol	Yohimbine
celecoxib	Etoricoxib	minoxidil	pregabalin	zolpidem
cetirizine	fluoxetine	mitragynine	prilocaine	methylone
chlorcyclizine	furosemide	modafinil	Propoxyphene	ethylone

Coverage vs Traditional Workplace test



Enhanced vs traditional workplace test on samples analysed to date



LGC Classification: NOT PROTECTIVELY MARKED



Population A

Specific study population – already undergoing ‘standard’ urine drug testing

1. Thousands of samples run
2. HRAM screening positive rate double the traditional type testing
3. Positives included controlled (classified) drugs not in routine testing panels – bufotenine
4. Non controlled substances which are suspected to have psychoactive outcomes not in routine testing – glaucine, 5F AKB-48.
5. Other substances (medication) available on prescription only not included in routine testing – zopiclone, gabapentin, pregabalin



Population B

Specific study population – already undergoing urine drug testing

1. Thousands of samples run
2. Current programme includes rudimentary NPS screening
3. HRAM screening positive rate over 18% on samples previously screened as negative
4. Anabolic steroid positive rate over 10%

Synthetic Cannabinoids detected in study



5F AB-PINACA	EAM-2201
5F AKB-48	FAM-2201
5F PB-22	MAM2201
5F UR-144	MDMB-CHMICA
AB-CHMINACA	STS-135
AM-2201	THJ-018
AM694	THJ-2201
Cumyl 5F PINACA	UR-144
BB-22	



LGC Classification: NOT PROTECTIVELY MARKED



Challenges

Challenges


- Standards
 - Often not available, purchase and purify
- Metabolites
 - Often none are published, In-vitro options
- Interpretation
 - Limited peer reviewed options
- Employment Policy
 - Often limited to 'impairment', controlled substances, or listed drugs
- Current
 - Keeping up to date, 'white powder' analysis



In Vitro Metabolism




- Knowledge of breakdown in the body is required – metabolism
- Many drugs only leave detectable concentrations of metabolites and not parent drug
- Metabolite information for many NPS are not available
- Able to mimic metabolism in a test tube using the drug and enzymes from liver cells
- Information is then used in testing process



LGC Classification: NOT PROTECTIVELY MARKED

MDMB CHMICA



- Current synthetic cannabinoid in UK
- Related to deaths across Europe and USA
- No metabolism information
- No reference material for metabolite


In vitro data for MDMB-CHMICA



\\INS-EXAC02\Data\...126Mar15\MDMBCHMICA
UK-LGC WORKPLACE TESTING

03/26/15 12:25:10

RT: 0.00 - 6.51




MS/MS spectra on MDMB-CHMICA hydroxy metabolite



MDMB-CHMICA_fixed_CE
UK-LGC WORKPLACE TESTING


04/01/15 10:11:34

RT: 0.00 - 6.50



NL: 1.09E9
TIC F: FTMS + c ESI
Full ms2
401.20@hcd35.00
[50.00-425.00] MS
MDMB-
CHMICA_fixed_CE

MDMB-CHMICA_fixed_CE #1553 RT: 4.04 AV: 1 NL: 7.57E8
T: FTMS + c ESI Full lock ms2 401.20@hcd35.00 [50.00-425.00]



Major metabolite in positive urine




H744645-1Cq_MDMB-CHMICA-mets
UK-NOMS-MDT

03/18/15 19:23:13


H744645-1C

RT: 0.00 - 6.50



NL: 3.10E7
TIC F: FTMS + c ESI
Full ms2
401.20@hcd40.00
[50.00-425.00] MS
H744645-1Cq_MDMB-
CHMICA-mets

H744645-1Cq_MDMB-CHMICA-mets #1556 RT: 4.03 AV: 1 SB: 1 3.23, 3.36 NL: 1.54E7
F: FTMS + c ESI Full ms2 401.20@hcd40.00 [50.00-425.00]





Questions?

- Can we use in vitro generated metabolites as reference materials where no other alternatives exist?